

Promote Accel™
CRT-D



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.



AUTOMATICALLY ADJUSTS TO CHANGING PATIENT NEEDS FOR OPTIMAL THERAPY DELIVERY.

The Promote Accel™ CRT-D from St. Jude Medical is designed to reduce life-threatening risks and simplify the implant, monitoring and follow-up processes for enhanced safety and efficiency. Featuring capture confirmation, biventricular trigger mode and AT/AF alerts, the Promote Accel CRT-D allows physicians to tailor therapy to changing patient needs to help ensure optimal therapy delivery and enable timely remote patient follow-up for safer, more efficient patient care.

Accelerating Patient Care

By providing ongoing patient monitoring and automatic threshold adjustments, the Promote Accel™ CRT-D is designed to offer faster, more efficient therapy delivery and management.

BiVCap™ Confirm

Automatic pacing adjustments meet changing patient needs.

- Automatically, periodically and independently measures and trends capture thresholds of both the right and left ventricles
- Minimizes risks associated with unpredictable changes in capture thresholds
- Provides reliable, automatic pacing adjustment



ACap™ Confirm

Automatic adjustments are designed to efficiently deliver optimal therapy.

- Utilizes an enhanced algorithm that increases the patient recommendation rate and does not require separate set-up tests
- Captures threshold trend data
- Provides reliable automatic pacing adjustment



A. BiVCap Confirm

LVCap™ Confirm and RVCap™ Confirm features ensure capture of the myocardium in response to pacing stimuli in the left and right ventricles.

B. ACap Confirm

ACap Confirm feature periodically completes a threshold search and automatically adjusts amplitudes to address patients' changing atrial thresholds.



BiV Trigger Mode

Triggered pacing ensures BiV therapy and maintains delivery of cardiac resynchronization therapy.

- Triggers pacing in response to a sensed event
- Ensures BiV therapy delivery during atrial arrhythmias

A. BiV Trigger

Triggered pacing with BiV Trigger Mode helps maintain a high percentage of BiV pacing by triggering pacing in response to a sensed event.



AT/AF Alerts

Alerts provide valuable insight into patient status.

- Provide an AT/AF burden alert when program threshold is exceeded
- Provide an AT/AF alert when continuous episode duration exceeds threshold
- Offer enhanced notification capabilities through vibratory patient notifier and Merlin™ Patient Care System

The **Promote Accel CRT-D** provides monitoring information and automatically adjusts output to meet changing patient needs to help ensure optimal therapy delivery and better patient outcomes.

St. Jude Medical is focused on reducing risk by continuously finding ways to put more control into the hands of those who save and enhance lives.

ATRIAL FIBRILLATION CARDIAC RHYTHM MANAGEMENT CARDIAC SURGERY CARDIOLOGY NEUROMODULATION

Global Headquarters

One Lillehei Plaza
St. Paul, Minnesota 55117
USA
+1 651 483 2000
+1 651 490 4310 Fax

Cardiac Rhythm

Management Division
15900 Valley View Court
Sylmar, California 91342
USA
+1 818 362 6822
+1 818 364 5814 Fax

St. Jude Medical

Sweden AB
Veddestavägen 19
175 84 Järfälla
Sweden
+46 8 474 4000
+46 8 760 9542 Fax

St. Jude Medical Japan Co., Ltd.

3-1-30, Minami-Aoyama
Minato-ku
Tokyo 107 0062
Japan
+81 3 3423 6450
+81 3 3402 5586 Fax

St. Jude Medical Europe, Inc.

The Corporate Village
Figuerras Building
Avenue Da Vinci laan, 11 Box F1
B-1935 Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical Brasil Ltda.

Rua Frei Caneca, 1380
7º ao 9º andares
01307-002 - São Paulo (SP)
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical (Hong Kong) Ltd.

Unit 2701-07 27/F, COSCO Tower
Grand Millennium Plaza
183 Queen's Road
Central, Hong Kong
+852 2996 7688
+852 2956 0622 Fax

sjm.com



Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country.

Unless otherwise noted, [™] indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are trademarks and service marks of St. Jude Medical, Inc. and its related companies.

©2009 St. Jude Medical, Inc. All Rights Reserved.

Item GMCRM390EN